

# veštačke proteze kuka



## synthetic hip prosthesis

## Particijalna proteza kuka Austin Moore / Partial hip prosthesis Austin Moore


| Kataloški broj<br>/ Catalog number | Prečnik glave [mm]<br>/ Head Diameter [mm] | Dužina vrata (mm)<br>/ Length of neck (mm) |
|------------------------------------|--------------------------------------------|--------------------------------------------|
| AM38                               | 38                                         |                                            |
| AM39                               | 39                                         |                                            |
| AM40                               | 40                                         |                                            |
| AM41                               | 41                                         |                                            |
| AM42                               | 42                                         |                                            |
| AM43                               | 43                                         |                                            |
| AM44                               | 44                                         |                                            |
| AM45                               | 45                                         |                                            |
| AM46                               | 46                                         |                                            |
| AM47                               | 47                                         |                                            |
| AM48                               | 48                                         |                                            |
| AM49                               | 49                                         |                                            |
| AM50                               | 50                                         |                                            |
| AM51                               | 51                                         |                                            |
| AM52                               | 52                                         |                                            |
| AM53                               | 53                                         |                                            |
| AM54                               | 54                                         |                                            |
| AM55                               | 55                                         |                                            |
| AM56                               | 56                                         |                                            |
| AM57                               | 57                                         |                                            |
| AM58                               | 58                                         |                                            |
| AM59                               | 59                                         |                                            |
| AM60                               | 60                                         |                                            |
| AM61                               | 61                                         |                                            |
| AM62                               | 62                                         |                                            |

34

Parcijalnim protezama kuka zamenjuje se glava i vrat butne kosti bez zamene acetabuluma. Ove proteze je u primenu uveo Amerikanac Austin Moore 1939. godine i još uvek su u upotrebi jer je operacija jednostavnija, kraće traje u odnosu na totalne endoproteze, i isplativije su. U principu se primenjuju kod osoba preko 65 godina.

Partial hip prosthesis replaces the head and femoral neck without replacing acetabulum. This prosthesis was introduced by American Austin Moore in 1939. It is still used because the operation is simpler, takes less time compared to total endoprosthesis and is cost effective. It is being used for people over 65 years of age.

Materijal: ISO 5832-1 (medicinski čelik)  
/ Material for implantation Stainless Steel ISO 5832/1




## Parcijalna proteza kuka (Austin Moore)

/Partial hip prosthesis "Austin Moore"

Uska  
/Narrow

| Kataloški broj<br>/Catalog number | Prečnik glave [mm]<br>/ Head Diameter [mm] | Dužina vrata (mm)<br>/Length of neck (mm) |
|-----------------------------------|--------------------------------------------|-------------------------------------------|
| AM38.U                            | 38                                         | 34                                        |
| AM39.U                            | 39                                         |                                           |
| AM40.U                            | 40                                         |                                           |
| AM41.U                            | 41                                         |                                           |
| AM42.U                            | 42                                         |                                           |
| AM43.U                            | 43                                         |                                           |
| AM44.U                            | 44                                         |                                           |
| AM45.U                            | 45                                         |                                           |
| AM46.U                            | 46                                         |                                           |
| AM47.U                            | 47                                         |                                           |
| AM48.U                            | 48                                         |                                           |
| AM49.U                            | 49                                         |                                           |
| AM50.U                            | 50                                         |                                           |

Materijal: ISO 5832-1 (medicinski čelik)  
/ Material for implantation Stainless Steel ISO 5832/1






## Parcijalna proteza kuka biartikularna /Partial biarticular hip prosthesis



### Telo parcijalne proteze kuka biartikularne / Body of partial biarticular hip prosthesis



Konus: 12/14 mm  
Taper: 12/14 mm

| Kataloški broj<br>/ Catalog number | Veličina tela - broj<br>/ Body size . - no. | Dužina tela [mm]<br>/ Body length [mm] | Dužina vrata tela [mm]<br>/ Neck length [mm] |
|------------------------------------|---------------------------------------------|----------------------------------------|----------------------------------------------|
| NA.TBP.1                           | 1                                           | 130.0                                  | 38.5                                         |
| NA.TBP.2                           | 2                                           | 137.5                                  |                                              |
| NA.TBP.3                           | 3                                           | 145.0                                  |                                              |
| NA.TBP.4                           | 4                                           | 152.5                                  |                                              |
| NA.TBP.5                           | 5                                           | 160.0                                  |                                              |
| NA.TBP.6                           | 6                                           | 167.5                                  |                                              |
| NA.TBP.7                           | 7                                           | 175.0                                  |                                              |
| NA.TBP.8                           | 8                                           | 182.5                                  |                                              |
| NA.TBP.9                           | 9                                           | 190.0                                  |                                              |

Materijal: ISO 5832-9 (medicinski čelik)  
/ Material for implantation Stainless Steel: ISO 5832-9



## Parcijalna proteza kuka sa izmenjivom glavom / Partial hip prosthesis with changeable head

### Izmenjiva glava biartikularne parcijalne proteze kuka /Changeable head


| Kataloški broj<br>/ Catalog number | Veličina<br>/ Size | Konus<br>/Taper |
|------------------------------------|--------------------|-----------------|
| IG 42                              | 42                 | 12/14           |
| IG 43                              | 43                 |                 |
| IG 44                              | 44                 |                 |
| IG 45                              | 45                 |                 |
| IG 46                              | 46                 |                 |
| IG 47                              | 47                 |                 |
| IG 48                              | 48                 |                 |
| IG 49                              | 49                 |                 |
| IG 50                              | 50                 |                 |
| IG 51                              | 51                 |                 |
| IG 52                              | 52                 |                 |
| IG 53                              | 53                 |                 |
| IG 54                              | 54                 |                 |
| IG 55                              | 55                 |                 |
| IG 56                              | 56                 |                 |
|                                    |                    |                 |
|                                    |                    |                 |



Materijal: ISO 5832-1 (medicinski čelik) I UHMWPE  
/ Material for implantation Stainless Steel: ISO 5832/1 and UHMWPE

glava parcijalne proteze kuka biartikularne-acetabulum  
/ head of partial biarticular hip prosthesis

| Kataloški broj<br>/ Catalog number | Prečnik [mm]<br>/ Diameter [mm] | Veličina prstena (mm)<br>/ Ring size (mm) |
|------------------------------------|---------------------------------|-------------------------------------------|
| NA.BP.40                           | 40                              | 31.5                                      |
| NA.BP.42                           | 42                              |                                           |
| NA.BP.44                           | 44                              |                                           |
| NA.BP.46                           | 46                              |                                           |
| NA.BP.48                           | 48                              |                                           |
| NA.BP.50                           | 50                              | 35.5                                      |
| NA.BP.52                           | 52                              |                                           |
| NA.BP.54                           | 54                              |                                           |
| NA.BP.56                           | 56                              |                                           |
| NA.BP.58                           | 58                              |                                           |



za glave/ for heads Ø32 mm; Ø28 mm

Umetak UHMWPE  
/ Insert: UHMWPEfemoralna glava  
/femoral head

| Kataloški broj<br>/ Catalog number | Veličina<br>/ Size                      | Prečnik (mm)<br>/ Diameter (mm) | Konus<br>/ Taper |
|------------------------------------|-----------------------------------------|---------------------------------|------------------|
| TCP 32GS                           | kratka/small                            | 32                              | 12/14            |
| TCP 32GM                           | srednja/medium                          |                                 |                  |
| TCP 32GL                           | duga/long                               |                                 |                  |
| TCP 32GXL                          | ekstra duga/extra long                  |                                 |                  |
| TCP 32GXXL                         | ekstra ekstra duga<br>/extra extra long |                                 |                  |

| Kataloški broj<br>/ Catalog number | Veličina<br>/ Size                      | Prečnik (mm)<br>/ Diameter (mm) | Konus<br>/ Taper |
|------------------------------------|-----------------------------------------|---------------------------------|------------------|
| TCP 28GS                           | kratka/small                            | 28                              | 12/14            |
| TCP 28GM                           | srednja/medium                          |                                 |                  |
| TCP 28GL                           | duga/long                               |                                 |                  |
| TCP 28GXL                          | ekstra duga/extra long                  |                                 |                  |
| TCP 28GXXL                         | ekstra ekstra duga<br>/extra extra long |                                 |                  |
| TCP 24GM                           |                                         | 24                              | 12/14            |



Materijal: ISO 5832-1

/ Material for implantation: ISO 5832/1