

VEŠTAČKE PROTEZE KUKA



SYNTHETIC HIP PROSTHESES

Particijalna proteza kuka "Austin Moore"

/ Partial hip prosthesis "Austin Moore"


Kataloški broj / Catalog number	Prečnik glave [mm] / Head Diameter [mm]	Dužina vrata (mm) / Length of neck (mm)
AM38	38	
AM39	39	
AM40	40	
AM41	41	
AM42	42	
AM43	43	
AM44	44	
AM45	45	
AM46	46	
AM47	47	
AM48	48	
AM49	49	
AM50	50	
AM51	51	
AM52	52	
AM53	53	
AM54	54	
AM55	55	
AM56	56	
AM57	57	
AM58	58	
AM59	59	
AM60	60	
AM61	61	
AM62	62	

34

Parcijalnim protezama kuka zamenjuje se glava i vrat butne kosti bez zamene acetabuluma. Ove proteze je u primenu uveo Amerikanac Austin Moore 1939. godine i još uvek su u upotrebi jer je operacija jednostavnija, kraće traje u odnosu na totalne endoproteze, i isplativije su. U principu se primenjuju kod osoba preko 65 godina.


Partial hip prosthesis replaces the head and femoral neck without replacing acetabulum. This prosthesis was introduced by American Austin Moore in 1939. It is still used because the operation is simpler, takes less time compared to total endoprosthesis and is cost effective. It is being used for people over 65 years of age.

Materijal: ISO 5832-1 (medicinski čelik)
/ Material for implantation Stainless Steel ISO 5832/1




Parcijalna proteza kuka "Austin Moore"
/Partial hip prosthesis "Austin Moore"**Uska**
/Narrow

Kataloški broj / Catalog number	Prečnik glave [mm] / Head Diameter [mm]	Dužina vrata (mm) / Length of neck (mm)
AM38.U	38	34
AM39.U	39	
AM40.U	40	
AM41.U	41	
AM42.U	42	
AM43.U	43	
AM44.U	44	
AM45.U	45	
AM46.U	46	
AM47.U	47	
AM48.U	48	
AM49.U	49	
AM50.U	50	

Materijal: ISO 5832-1 (medicinski čelik)
/ Material for implantation Stainless Steel ISO 5832/1

Totalna cementna proteza kuka /Cemented total hip prosthesis

telo totalne cementne proteze kuka /Body of cemented total hip prosthesis



telo totalne cementne proteze kuka body of cemented total hip prosthesis

Kataloški broj / Catalog number	Veličina tela - broj / Body size - number	Dužina tela (mm) / Body length (mm)	Debljina (mm) / Thickness (mm)	Konus /Taper
TCPT1	1	141	10.5	12/14 ili/or 14/16
TCPT2	2	145	11.5	12/14 ili/or 14/16
TCPT3	3	157	12	12/14 ili/or 14/16
TCPT4	4	168	13	12/14 ili/or 14/16
TCPT5	5	178	14	12/14 ili/or 14/16

Materijal: ISO 5832-9 (medicinski čelik)
/ Material for implantation Stainless Steel: ISO 5832/9

Totalna cementna proteza kuka /Cemented total hip prosthesis

glava totalne cementne proteze kuka /Head of cemented total hip prosthesis

Kataloški broj / Catalog number	Veličina / Size	Prečnik (mm) / Diameter (mm)	Konus /Taper
TCP 28GS	kratka/small	28	12/14
TCP 28GM	srednja/medium		
TCP 28GL	duga/long		
TCP 28GXL	ekstra duga/extra long		
TCP 28GXXL	ekstra ekstra duga /extra extra long		



Materijal: ISO 5832-9 (medicinski čelik)
/ Material for implantation Stainless Steel: ISO 5832/9

Kataloški broj / Catalog number	Veličina / Size	Prečnik (mm) / Diameter (mm)	Konus /Taper
TCP 32GS	kratka/small	32	14/16
TCP 32GM	srednja/medium		
TCP 32GL	duga/long		
TCP 32GXL	ekstra duga/extra long		
TCP 32GXXL	ekstra ekstra duga /extra extra long		



Materijal: ISO 5832-9 (medicinski čelik)
/ Material for implantation Stainless Steel: ISO 5832/9

acetabulum totalne cementne proteze kuka /acetabulum of cemented total hip prosthesis

za glave/ for heads Ø32 mm; Ø28 mm

Kataloški broj / Catalog number	Veličina / Size
TPAC 46	46
TPAC 48	48
TPAC 50	50
TPAC 52	52
TPAC 54	54
TPAC 56	56
TPAC 58	58
TPAC 60	60
TPAC 62	62




Materijal: ISO 5834-2 (medicinska plastika)
/ Material for implantation medical plastic: ISO 5834/2

NAPOMENA: Pored dimenzija navedenih u tabelama, mogu se izrađivati i po posebnim zahtevima
NOTE: Additional dimensions can be made upon request

Parcijalna proteza kuka sa izmenjivom glavom / Partial hip prosthesis with changeable head

Telo parcijalne proteze kuka sa izmenjivom glavom / Body of the partial hip prosthesis with changeable head



Kataloški broj / Catalog number	Veličina tela - broj / Body size - number	Dužina tela (mm) / Body length (mm)	Debljina (mm) / Thickness (mm)	Konus /Taper
TCPT1	1	141	10.5	12/14 ili/or 14/16
TCPT2	2	145	11.5	12/14 ili/or 14/16
TCPT3	3	157	12.0	12/14 ili/or 14/16
TCPT4	4	168	13.0	12/14 ili/or 14/16
TCPT5	5	178	14.0	12/14 ili/or 14/16

Materijal: ISO 5832-9 (medicinski čelik)
/ Material for implantation Stainless Steel: ISO 5832/9

Parcijalna proteza kuka sa izmenjivom glavom / Partial hip prosthesis with changeable head

Izmenjiva glava parcijalne proteze kuka /Changeable head

Kataloški broj / Catalog number	Veličina / Size	Konus /Taper
IG 38	38	12/14
IG 39	39	
IG 40	40	
IG 41	41	
IG 42	42	
IG 43	43	
IG 44	44	
IG 45	45	
IG 46	46	
IG 47	47	
IG 48	48	
IG 49	49	
IG 50	50	
IG 51	51	
IG 52	52	
IG 53	53	
IG 54	54	
IG 55	55	
IG 56	56	
IG 57	57	
IG 58	58	
IG 59	59	
IG 60	60	
IG 61	61	
IG 62	62	



Materijal: ISO 5832-9 (medicinski čelik)
/ Material for implantation Stainless Steel: ISO 5832/9

Parcijalna proteza kuka biartikularna /Partial biarticular hip prosthesis



Telo parcijalne proteze kuka biartikularne / Body of partial biarticular hip prosthesis



Kataloški broj / Catalog number	Veličina tela - broj / Body size . - no.	Dužina tela [mm] / Body length [mm]	Dužina vrata tela [mm] / Neck length [mm]
NA.TBP.1	1	130.0	38.5
NA.TBP.2	2	137.5	
NA.TBP.3	3	145.0	
NA.TBP.4	4	152.5	
NA.TBP.5	5	160.0	
NA.TBP.6	6	167.5	
NA.TBP.7	7	175.0	
NA.TBP.8	8	182.5	
NA.TBP.9	9	190.0	

Materijal: ISO ISO 5832-9 (medicinski čelik)
/ Material for implantation Stainless Steel: ISO 5832-9

glava parcijalne proteze kuka biartikularne
/ head of partial biarticular hip prosthesissklop
/ assembledbiartikularna glava
/ biarticular head

Kataloški broj / Catalog number	Prečnik [mm] / Diameter [mm]	Veličina prstena (mm) / Ring size (mm)	Prečnik glavice (mm) / Head diameter (mm)
NA.BP.42	42	31.5	24
NA.BP.44	44		
NA.BP.46	46		
NA.BP.48	48		
NA.BP.50	50	35.5	28
NA.BP.52	52		
NA.BP.54	54		
NA.BP.56	56		

Materijal (umetak) ISO 5834-2
/ Material for implantation(insert): ISO 5834/2Plajt: ISO 5832-1 (medicinski čelik)
/ Insert: ISO 5832-1 (Stainless Steel)femoralna glava
/ femoral head

Kataloški broj / Catalog number	Veličina / Size	Prečnik (mm) / Diameter (mm)	Konus / Taper
TCP 28GS	kratka/small	28	12/14
TCP 28GM	srednja/medium		
TCP 28GL	duga/long		
TCP 28GXL	ekstra duga/extra long		
TCP 28GXXL	ekstra ekstra duga /extra extra long	24	12/14
TCP 24GM			

Materijal : ISO 5832-9
/ Material for implantation: ISO 5832/9Materijal: ISO 5834-2
/ Material for implantation: ISO 5834/2prsten
ring